

# CBC-CAL PLUS

## HEMATOLOGY CALIBRATOR

### CALIBRATOR


**LOT** PLUS0625  
 **2025-07-05**

Table of Assigned Values <sup>1</sup>	WBC 10 <sup>3</sup> /μL	RBC 10 <sup>6</sup> /μL	HGB g/dL	HCT %	MCV fL	PLT 10 <sup>3</sup> /μL	MPV fL
<b>DxH™ Series Diluent</b> Beckman Coulter® UniCel® DxH™ 900/800/690T/600	9.0	4.42	13.2	38.9	88.0	232	7.7
<b>Tolerance Limits (±)</b>	0.2	0.05	0.1	0.8	1.0	10	0.5

Record Results and Calculate Mean Value	WBC 10 <sup>3</sup> /μL	RBC 10 <sup>6</sup> /μL	HGB g/dL	HCT %	MCV fL	PLT 10 <sup>3</sup> /μL	MPV fL
<b>Sample <sup>2</sup></b> 1							
2							
3							
4							
5							
6							
7							
8							
9							
10							
Mean Value							

Determine Adjustment Factor <sup>3</sup>	WBC 10 <sup>3</sup> /μL	RBC 10 <sup>6</sup> /μL	HGB g/dL	HCT %	MCV fL	PLT 10 <sup>3</sup> /μL	MPV fL
<u>Assigned Value</u>	_____	_____	_____	_____	_____	_____	_____
Mean Value	_____	_____	_____	_____	_____	_____	_____
Calibration Adjustment Factor (round off to 3 decimal places)							

Record Verification Results	WBC 10 <sup>3</sup> /μL	RBC 10 <sup>6</sup> /μL	HGB g/dL	HCT %	MCV fL	PLT 10 <sup>3</sup> /μL	MPV fL
<b>Sample <sup>2</sup></b> 1							
2							
3							
4							
5							
6							
7							
8							
9							
10							
Mean Value							

Before using, refer to the instruction sheet for mixing directions.

For technical assistance in the USA and Canada call Technical Service at (800) 523-3395.

All brands and products are trademarks or registered trademarks of their respective companies.

#### NOTES:

If analyzer provides both a corrected and uncorrected WBC count, use uncorrected WBC count for calibration purposes.

Disregard the HCT or MCV value if it is not used to calibrate your system.

Analyze an appropriate number of samples for your system. A minimum of six is recommended.

Determine a calibration adjustment factor if required for your system.

**R&D Systems, Inc.**  
 614 McKinley Place NE  
 Minneapolis, MN 55413 USA  
 AS018-021 Rev. 04/23