

CBC-CAL PLUS

For PROKAN HEMATOLOGY CALIBRATOR

ASSAY VALUES AND EXPECTED RANGES

LOT**PLUS0421****2021-05-05**

Instrument	Parameter	Target Value	Acceptable Range
PROKAN PE-6000	WBC $\times 10^3/\mu\text{L}$	9.0	± 0.2
	RBC $\times 10^6/\mu\text{L}$	4.65	± 0.08
	HGB g/dL	14.1	± 0.2
	MCV fL	86.0	± 2.0
	PLT $\times 10^3/\mu\text{L}$	251	± 15
PROKAN PE-6100	WBC $\times 10^3/\mu\text{L}$	9.0	± 0.2
	RBC $\times 10^6/\mu\text{L}$	4.65	± 0.08
	HGB g/dL	14.1	± 2
	MCV fL	86.0	± 2.0
	PLT $\times 10^3/\mu\text{L}$	251	± 15
PROKAN PE-6800	WBC $\times 10^3/\mu\text{L}$	9.1	± 0.2
	RBC $\times 10^6/\mu\text{L}$	4.70	± 0.08
	HGB g/dL	14.0	± 2
	MCV fL	86.7	± 2.0
	PLT $\times 10^3/\mu\text{L}$	253	± 15

WBC Differential alarms may occur with calibrator material. This will not affect the validity of results.

Before using, refer to the instruction sheet for mixing directions.

Shenzhen Prokan Electronics Inc.

No.1002, Building 4, Fantasia MIC Plaza, Pengji Pioneering Park,
the west of Nanhai Road, Nanshan District, Shenzhen, China
TEL: +86-755-26955166 FAX: +86-755-26747910

CBC-CAL PLUS

For PROKAN
HEMATOLOGY CALIBRATOR

ASSAY VALUES AND EXPECTED RANGES

LOT PLUS0421
 2021-05-05

Instrument	Parameter	Target Value	Acceptable Range
PROKAN PE-7100	WBC $\times 10^3/\mu\text{L}$	8.9	± 0.2
	RBC $\times 10^6/\mu\text{L}$	4.50	± 0.08
	HGB g/dL	13.7	± 0.2
	MCV fL	93.4	± 2.0
	PLT $\times 10^3/\mu\text{L}$	245	± 15

WBC Differential alarms may occur with calibrator material. This will not affect the validity of results.

Before using, refer to the instruction sheet for mixing directions.

Shenzhen Prokan Electronics Inc.

No.1002, Building 4, Fantasia MIC Plaza, Pengji Pioneering Park,
the west of Nanhai Road, Nanshan District, Shenzhen, China

TEL: +86-755-26955166 FAX: +86-755-26747910