

CD-Cal

HEMATOLOGY CALIBRATOR

CALIBRATOR

LOT

3KC0525



2025-06-05

Method	Parameter	Calibration Values	Acceptable Range
Cell-Dyn Sapphire™	WBC 10 ⁹ /L	8.80	8.60 - 9.00
	RBCi 10 ¹² /L	4.45	4.40 - 4.50
	RBCo 10 ¹² /L	4.45	4.40 - 4.50
	HGB g/dL	13.2	13.0 - 13.4
	MCV fL	79.5	78.5 - 80.5
	PLTo 10 ⁹ /L	215	205 - 225
	PLTi 10 ⁹ /L	237	227 - 247
	MPV fL	8.00	7.50 - 8.50
Cell-Dyn 3200 Software Version 1.8 or higher	WBC (WOC) 10 ⁹ /L	8.9	8.7 - 9.1
	NOC 10 ⁹ /L	8.9	8.7 - 9.1
	RBC 10 ¹² /L	4.38	4.32 - 4.44
	HGB g/dL	13.0	12.8 - 13.2
	MCV fL	76.0	75.0 - 77.0
	PLT 10 ⁹ /L	221	211 - 231
Cell-Dyn Ruby™	WBC (WOC) 10 ⁹ /L	8.9	8.7 - 9.1
	NOC 10 ⁹ /L	8.9	8.7 - 9.1
	RBC 10 ¹² /L	4.38	4.32 - 4.44
	HGB g/dL	13.0	12.8 - 13.2
	MCV fL	76.0	75.0 - 77.0
	PLT 10 ⁹ /L	221	211 - 231

Before using, refer to the instruction sheet for mixing directions.

For technical assistance in the USA and Canada call Technical Service at (800) 523-3395.

All brands and products are trademarks or registered trademarks of their respective companies.